


Order of Draw

	<ul style="list-style-type: none"> • Blood Culture • Clear Plastic Tube • Blue Verify Now Tube • Light Blue Tube • Royal Blue Red Label • Red (Clot Activator) • Gold (Clot Activator) • Light or Dark Green Lithium Heparin Tube • Dark Green Sodium Heparin Tube • Pink EDTA tube • Lavender EDTA tube • Royal Blue Lavender Label • Tan • Gray • Yellow ACD (AorB) Tube • VBG syringe
--	--