


# Microbiology Collection Containers

## BACTERIAL, FUNGAL (YEAST) COLLECTION KIT



**USE FOR:**  
 Routine Culture (Wounds),  
 MRSA PCR,  
 Group B strep PCR,  
 Throat Strep Screen,  
 Fungus (Yeast),  
 Vaginitis Exam  
*Trichomonas* antigen


## VIRAL COLLECTION KIT



**USE FOR:**  
 ED/Urgent Care Influenza,  
 RSV, PCR  
 Inpatient/Outpatient Influenza,  
 RSV, PCR  
 Respiratory Viral Panel  
 Expanded, PCR  
 Herpes Simplex Virus, PCR  
 Varicella Zoster Virus, PCR

## STERILE CONTAINER

Do not use if the blue seal is broken



**USE FOR:**  
 Urine Culture  
 Tissue Culture  
 C. difficile PCR  
 Fecal Lactoferrin for WBCs  
 Fungal Culture  
 H. pylori Stool Antigen,  
 Mycobacterial Culture  
 Respiratory Culture


## ANAEROBIC TRANSPORT VIAL

Do not use if vial has been removed from sterile packaging



**USE FOR:**  
 Anaerobic Cultures for  
 Bone  
 Fluid  
 Fungal  
 Mycobacterial  
 Tissue

## BACTERIAL COLLECTION KIT (PARA PAK CAREY BLAIR)



**USE FOR:**  
 Fecal Bacterial Pathogen PCR  
 Feces Culture unusual Pathogens

## TOTAL FIX VIAL



**USE FOR:**  
 Ova and Parasite  
 Giardia and Cryptosporidium Exam

Microbiology Instructions for collection of Aptima vials (GC and Chlamydia) and Virology collection are also available. Contact Laboratory Customer Service for information. or for information about specific test collection not mentioned here

**GO TO** our test catalog

<https://UVMLabs.TestCatalog.org/> or call 847-5121