


SPECTRUM HEALTH





Quality In Quality Out

**Labeling
Microtainer
Specimens**



...or how to label those itty bitty tubes





Understanding the
Seven Part Cerner Label
will help

What is the 7 part Cerner label?



Making more sense?



Main Label **1**

Aliquot Labels **2**

Aliquot Labels **3**

Aliquot Labels **4**

Aliquot Labels **5**


Barcode Labels **6**

Barcode Labels **7**

Ta Da !

001646136 Blodgett		DSCH	
CRNTEST, BLODGETT		53 Y F AUG-15-52	
		NOV-01-05 1200	
BC0016461367			
			
ACC#	05-305-03229	Coll List: 00	
BLUE (A)	RT/Ro	3mL LAV (B)	RT/Ro
Protime		CBC w/ AutoDiff	
DR. Alter, David N		DR. Alter, David N	
Deliver to: COAG		Deliver to: HEME	
001646136	001646136	001646136	001646136
CRNTEST, BLO	CRNTEST, BLO	CRNTEST, BLO	CRNTEST, BLO
DSCH	DSCH	DSCH	DSCH
NOV-01-05 1200	NOV-01-05 1200	NOV-01-05 1200	NOV-01-05 1200
05-305-03229A	05-305-03229A	05-305-03229B	05-305-03229B
Protime	Protime	CBC w/ Aut	CBC w/ Aut
			
001646136 AUG-15-52 FC0AG	001646136 AUG-15-52 FH0ME		
CRNTEST, BLODGETT	CRNTEST, BLODGETT		
DSCH	DSCH		
NOV-01-05 1200	NOV-01-05 1200		
05-305-03229A	05-305-03229B		
Protime	CBC w/ AutoDiff		

But what goes where?




Barcode vs. Aliquot Labels



HEMATOLOGY / BLOOD BANK

BARCODE LABELS

- *LAVENDER* color
- *VERTICAL* placement
- *BELOW* cap
- *NOT* wrapped or wrinkled




Correct

Incorrect

COAGULATION

BARCODE LABELS

- *BLUE* color
- *THREE* sizes
- *VERTICAL* labeling on all tubes



Correct



CHEMISTRY / IMMUNOLOGY

ALIUQUOT LABELS

- Green, yellow, red
- Place *ALIUQUOT* labels on specimen
- Send barcode labels with specimen.



Correct

Incorrect


Correct labeling is important



because.....

.....it provides

optimum
turnaround time
for lab results





Quality In = Quality Out

SPECTRUM HEALTH

