



LAB ALERT

RL Ryan H. Livengood, M.D., Medical Director
Department of Pathology and Laboratory Services

Facility: Berkeley Medical Center and Jefferson Medical Center

Effective: Immediately

New swab for collection of throat sample for Strep A molecular testing:



Effective immediately throat swabs for the Strep A molecular test should be collected using a Puritan foam applicator tipped swab. The particular swab is necessary to comply with manufacturer specifications.

*If you have any questions call BMC at 304-264-1212, option #7 or JMC at 304-728-1759