

Office/Facility name _____ Date _____

Address _____ Ordered by _____ Phone _____

PHLEBOTOMY			REQUISITIONS, FORMS AND PATIENT EDUCATION		
QTY			QTY		
	Tube gold 5 mL (SST)	(order as ea)		Laboratory supply requisition	(order as ea)
	Tube red 6 mL (plain)	(order as ea)		General lab requisition form #443 or #435	(100/pkg)
	Tube lavender 3 mL (EDTA)	(order as ea)		Anatomical pathology requisition form #447	(100/pkg)
	Tube light blue 2.7 mL (3.2% citrate)	(order as ea)		Advanced technology requisition form #480	(100/pkg)
	Tube light blue 4.5 mL (3.2% citrate)	(order as ea)		Microbiology requisition form # 490	(100/pkg)
	Tube dark green 6 mL (sod heparin)	(order as ea)		Cystic fibrosis requisition form #450	(100/pkg)
	Tube light green 4.5 mL (lith heparin PST)	(order as ea)		Allergen test requisition form #495	(100/pkg)
	Tube royal blue 6 mL (trace element)	(order as ea)		Toxicology requisition form #485	(100/pkg)
	Tube royal blue 6 mL mayo (metal free EDTA)	(order as ea)		Stain worksheet	(order as ea)
	Tube royal blue 6 mL mayo (no additive aluminum)	(order as ea)		Additional test requisition	(order as ea)
	Tube tan 3 mL (EDTA)	(order as ea)		HIV consent form (English)	(order as ea)
	Tube pink 6 mL (EDTA)	(order as ea)		HIV consent form (Spanish)	(order as ea)
	Tube gray 6 mL (sod fluoride)	(order as ea)		Cystic fibrosis brochure	(order as ea)
	Tube yellow 6 mL ACD (sol A)	(order as ea)		Quad PT education booklet (English)	(order as ea)
	Tube yellow 6 mL ACD (sol B)	(order as ea)		Quad PT education booklet (Spanish)	(order as ea)
	Tube microtainer green (lith heparin)	(25/bag)		Informed consent for genetic testing booklet	(order as ea)
	Tube microtainer gold (SST)	(25/bag)		Patient home collection instructions (English)	(order as ea)
	Tube microtainer red (plain)	(25/bag)		Patient home collection instructions (Spanish)	(order as ea)
	Tube microtainer lavender (EDTA)	(25/bag)		eSHare lab order slips (green)	(order as ea)
	Tube QuantiFeron	(order as ea)		Epic lab order slips (orange)	(order as ea)
	Tube urine myoglobin (mayo)	(order as ea)	QTY	CYTOLOGY	
	Vacutainer holders	(order as ea)		Pap test thin prep vial w/ cervical broom	(25/pkg)
	Multi draw vacutainer needle 21 - 1¼"	(48/box)		Pap test thin prep vial w/ cervical brush & spatula	(25/pkg)
	Multi draw vacutainer needle 22 - 1¼"	(48/box)		Fine needle aspirate (FNA) kit	(order as ea)
	Clear aliquot tube	(25/bag)		Cytolyt for non GYN small mouth	(order as ea)
	Brown/Amber aliquot tube	(25/bag)		Cytolyt for non GYN large mouth	(order as ea)
	Mayo metal free aliquot tube	(25/bag)		Green top slide fixative	(order as ea)
	Mayo clear aliquot tube	(100/bag)		UroVysion kits	(order as ea)
	Blood transfer device	(order as ea)	QTY	GLUCOLA	
	Blood culture device	(order as ea)		Glucola orange (50 gram)	(6/pkg)
QTY	TOXICOLOGY			Glucola orange (100 gram)	(6/pkg)
	Blue top drug screen cups	(order as ea)		Glucola varied flavors (100 gram)	(6/pkg)
	Drug screen cup seals (numbered)	(order as ea)	QTY	SPECIMEN BAGS	
	Drug screen box mailers	(order as ea)		White specimen bag (ambient)	(100/pkg)
	Drug screen seals (unnumbered)	(roll)		Green specimen bag (refrigerated)	(100/pkg)
	Drug screen box mailer address labels	(roll)		Blue specimen bag (frozen)	(100/pkg)
QTY	MICROBIOLOGY			Red specimen bag (stat)	(100/pkg)
	Standard tip Eswab	(order as ea)		Purple specimen bag (priority)	(100/pkg)
	Nasopharyngeal Eswab	(order as ea)		Biohazard specimen bag (large)	(order as ea)
	Blood culture bottle aerobic	(order as ea)		Biohazard specimen bag (extra large)	(order as ea)
	Blood culture bottle anaerobic	(order as ea)	QTY	MISCELLANEOUS	
	Blood culture prep kit	(order as ea)		Breath test kit	(order as ea)
	Hologic aptima swab kits vaginal #301162	(50/box)		Specimen log book	(order as ea)
	Hologic aptima urine kits M/F #301040	(order as ea)		Pediatric lead kit	(order as ea)
	Pinworm paddle	(order as ea)		72 hour stool container	(2/pkg)
	C & S, O & P stool kit	(order as ea)		Fetal fibronectin collection kit	(order as ea)
	STD multi-collection kit	(order as ea)		Panorama prenatal screen	(order as ea)
	UTM viral transport tube	(order as ea)	QTY	HISTOLOGY	
QTY	URINE COLLECTION			Container with formalin 20 mL	(order as ea)
	Sterile orange top urine container	(order as ea)		Container with formalin 40 mL	(order as ea)
	24 hour urine container	(order as ea)		Container with formalin 60 mL	(order as ea)
	Urine hat	(order as ea)		Container with formalin 180 mL	(order as ea)
	White 24 hour urine container (do not refrigerate)	(order as ea)		Container with formalin 480 mL	(order as ea)
	Urine C & S kit	(order as ea)		Container dry 180 mL	(order as ea)
	Pediatric urine collection kit (PUCK)	(order as ea)		Container dry 32 oz	(order as ea)
	Acetic acid (15 mL)	(6/pkg)		Container dry 83 oz	(order as ea)
	Acetic acid (25 mL)	(6/pkg)		Container dry 1 gallon	(order as ea)
	Boric acid 4-2.5 gram tablets	(order as ea)		Container dry 1 gallon/lids	(order as ea)
	Sodium carbonate 5 grams	(order as ea)	QTY	OTHER	
	Sodium hydroxide 10 mL	(order as ea)			
	Urine yellow vacutainer 10 mL	(order as ea)			
	Urine transfer straw	(order as ea)			
QTY	LABELS				
	7 part 1" core Cerner labels	(roll)			
	7 part 3" core Cerner labels	(roll)			
	Manual specimen labels	(roll)			
	eSHare printer labels	(roll)			

SEND ORIGINAL WITH YOUR COURIER OR FAX TO 616.267.2085
ALLOW 5 DAYS FOR DELIVERY OF SUPPLIES
Spectrum Health Laboratory Services Call-Center 616.774.7721