

BELLIN HEALTH OUTREACH LABORATORY
SUPPLY WORKSHEET

FACILITY: _____ DATE: _____ FILLED BY: _____

SUPPLY	AMOUNT	UNITS (CIRCLE)	INV #	SUPPLY	AMOUNT	UNITS (CIRCLE)	INV #
COLLECTION TUBES				FORMS			
LAVENDER 4.0 ML	BOX EACH		001 / 068	BLUE OUTREACH	PACK		031
BLUE 1.8 ML	BOX EACH		003 / 070	RED MC OUTREACH	PACK		033
GREEN PST 3.0 ML	BOX EACH		004 / 071	HISTOLOGY/TISSUE	PACK		034
GOLD SST 4.0 ML	BOX EACH		005 / 072	PAP DUPLICATE	PACK (100)		035
RED PLAN 6.0 ML	BOX EACH		006 / 073	PAP SINGLE SHEET/PAD	PAD (100)		089
PINK BLD BANK EDTA	BOX EACH		100 / 092	CYTOLOGY NON/GYN	PACK (100)		090
MICROTAINER LAV	BAG(50) EACH		012 / 013				
MICROTAINER BRN	BAG(50) EACH		014 / 015	CONTAINERS			
URINE CULTURE	BOX (50) EACH		016 / 093	FORMALIN 1 OZ	BOX (12)		036
URINALYSIS	BOX (50)		099 / 094	FORMALIN 2 OZ	BOX (12)		037
				URINE CUPS(STERILE)	BAG (100)	EACH	038 / 095
COLLECTION SUPPLIES				URINE MIDSTREAM KITS	BOX (50)	EACH	039 / 052
BD NEEDLE HOLDER	BAG (250) EACH		083 / 017	SPECIMEN TRANSPORT TUBES	BAG (100)	EACH	040 / 053
BD TRANSFER DEVICE	EACH		025	24 HR URINE CONTAINER	EACH		041
NEEDLE 22 X 1.25	BOX OF 100		019	SPEC TRANSPORT BAGS	PACK (100)		042
NEEDLE 21 X 1.25	BOX OF 100		021	TRANSFER PIPETTES	BOX (500)		043
BLD CULT BOTTLE -ANA	EACH		029				
BLD CULT BOTTLE -AER	EACH		030	MICROBIOLOGY			
BLD CULT BOTTLE - PED	EACH		032	AEROBE CULTURETTES	BAG(50) EACH		044 / 045
BUTTERFLY 21g	BOX OF 50		024	M5 MEDIA	BOX(72) EACH		046 / 047
BUTTERFLY 23g	BOX OF 50		023	DACRON/RAYON SWABS	BOX (100) EACH		048 / 049
BLOOD CULT PREPS	EACH		075	NASOPHARYNGEAL SWABS	EACH		050
QUANTIFERON KITS	EACH		103	URETHRAL SWABS FOR PCR	EACH		051
				STOOL C & S KITS	EACH		054
MISCELLANEOUS				STOOL O & P VALS	EACH		066
PAP KITS (1 SLIDE)/(2 SLIDE)	EACH		057 / 091	PINWORM PADDLES	EACH		056
TONER CARTRIDGE (SEE NOTE)	EACH		060	SPUTUM BULLETS	EACH		076
PAP KITS THIN PREP	EACH		085	BENZ-CHLORIDE WIPES	BOX		055
FROSTED SLIDES	BOX (50)		058	ACETIC ACID ALIQUOTS	EACH		064
SPECIMEN LABELS(SEE NOTE)	ROLL		059				
PRINTER PAPER (SEE NOTE)	REAM		062				

NOTES:

ALL SUPPLIES TO BE USED FOR BELLIN PERFORMED TESTING ONLY

Note: supplied for Bellin remote printer users only
Customer Service: 920-433-3652 Toll Free: 1-800-236-1639

Revised 5/16/11
10-6543draft