

Supplies

The following supplies are available through the Holland Hospital Laboratory Services packaging area. Please use our special supply requisitions (sent monthly via fax or courier) when ordering or you may call 616-355-3801.

For Blood Collection

Blood transfer device

Catecholamine tube, EDTA-sodium metabisulfite solution—

 10 mL for catecholamine fractionation, plasma, free

Finger puncture device

Gloves

 Extra small

 Small

 Medium

 Large

 X-Large

Glove sizes (Sterile)

 6 1/2

 7

 7 1/2

 8

 8 1/2

Heel puncture device

Needles

 21 gauge

 22 gauge

 21 gauge Butterfly

 23 gauge Butterfly

 25 gauge Butterfly

Needle holder

Serotonin collection tube

Syringes

 5 mL

 10 mL

 12 mL

VACUTAINERS®

 Gold (serum gel)—5 mL

 Green-top (sodium heparin)—6 mL

 Grey-top (potassium oxalate/sodium fluoride)—6 mL

 Lavender-top (EDTA)—2 mL and 4 mL

 Lavender-top MICROAINER® (EDTA)—0.5 mL

 Light blue-top (sodium citrate)—2.7 mL

 Mint-top (lithium heparin)—4.5 mL

 Pink-top (EDTA)—6 mL

 Red-top (plain)—6 mL

 Royal blue-top metal-free Monoject® (EDTA additive)—
 7 mL

 Royal blue-top metal-free Monoject® (no additive)—
 7 mL

MICROAINER® (serum gel tube)—0.5 mL

Yellow-top (ACD [solution B])—6 mL

MICROAINER® (lithium heparin)—0.5 mL

For Chromosome and/or Tissue Collection

Cytology fixative

Hank's balanced salt solution for chromosome analyses

Metal-free specimen vial (blue label) for tissue metal
 collections

Vial, 10% formalin

Zeus Transport Media (Michelle's Solution)

For Microbiology Testing

Anaerobic transport vial—liquid specimens—5 mL

Aerobic blood culture bottles (blue)—Adult

Anaerobic blood culture bottles (red)—Adult

Cobas® PCR Female Swab Sample Kit

Cobas® PCR Urine Sample Kit

ESwab™ collection kit

 for Aerobic and Anaerobic cultures

ESwab™ Nasopharyngeal collection kit

Formalin-Meridian 10% buffered neutral

Para-Pak® EcoFix® for Giardia and O & P

Pediatric blood culture bottle

Pinworm Paddles

Stool culture transport media

ThinPrep® Collection Kit

Tissue specimens—Sterile cup

UTM Viral Transport Media for Chlamydia Culture

Forms/Brochures

ABN Forms

Anatomic Pathology Requisition

Clinical Requisition

Cystic Fibrosis Brochures

1-hour Gestational Glucose Instructions

Lead Requisition

PSC Maps

Quad/Prenatal Serum Screening Requisition

Request for Supplies

Miscellaneous

Bags, specimen
Cleansing pads (Benzalkonium Chloride Towelettes)
FISH Kits
Flow media
FNA Kits
10% formalin bottles (120 mL, 60 mL, 20 mL)
Glucose tolerance beverage
Occult Blood Collection Kit
12-part Prostate Biopsy Kit
Pediatric urine collection bags
Semen Analysis Kits
Serum vial (6 mL), white cap
Stool container, large (24-hour)
Urine container, (24-hour graduated)
Urine cups, sterile
Urine collection hats

Holland Hospital Laboratory Services Test Catalog