


# Which Royal Blue tube should I use?

**BD vacutainer, trace element, EDTA, plastic tube, blue label with lavender strip**

This tube is for:

Arsenic, Cadmium, Heavy Metals, Magnesium RBC, Mercury



*Special order tube (order through Courier Supply Requisition)*

**BD vacutainer, trace element, no additive, serum plastic tube with red vertical strip on label**

This tube is for:

Chromium, Copper, Selenium, Cobalt, Iodine, Manganese, Zinc



*Available through SH Supply Chain Lawson #457005*

**Trace Element, No Additive, Aluminum Only**

This tube is for:

Aluminum

Manufacturers change frequently.  
**Please check tube label and lab catalog for correct tube.**  
If you are still unsure, contact the Referral's Department.



*Special order tube (order through Courier Supply Requisition)*