







Vaginal Specimen Collection Guide 47.05.T02

Test	Specimen/Media	Lawson Number
<p><u>Wet Prep</u>- microscopically looks for clue cells, yeast, and trichomonas.</p>	<p>Culturette swab</p>  <p>Black Top – Saline</p> 	<p>6117712</p> <p>6193366</p>
<p><u>Trichomonas PCR</u></p> <p>*Performed at South University Molecular Lab.</p>	<p>GeneXpert Swab</p> 	<p>6272690</p>
<p><u>Gonorrhea/Chlamydia</u></p> <p>*Performed at South University Molecular Lab.</p>	<p>Cobas collection kit</p> 	<p>6292565</p>
<p><u>Herpes</u></p> <p>*Performed at South University Molecular Lab.</p>	<p>Universal Transport Media</p> 	<p>6260951</p>
<p><u>Vaginitis Probe-PCR</u> testing for clue cells, yeast, and trichomonas. *Performed at Sanford Sioux Falls Reference Lab.</p>	<p>Affirm kit</p> 	<p>6185788</p> <p>*limited stock available from lab.</p>